

Sterile Collection Container



Ordered test(s):

C. Difficil
Chlam/GC DNA
Fecal leukocytes
Fungus Culture
(Nail, Skin, Hair)

Stool Culture
Sputum Culture
Synovial Fluid
Trichomonas
Urine culture

Stool Specimens



Ordered Test(s):
Ova and Parasites



Ordered Test(s):
Occult Blood

SPS Tube



Ordered Test(s):
Blood Culture

Clear Scotch Tape



Ordered Test(s):
Pinworm

Culture Swab



Ordered test(s):

GC Screen
Gram Stain
Rapid Strep A (2 swabs) -Throat
Strep Group B Screen-Vag
Throat, Vaginal, Ear, Eye,
Nasal, Rectal, Wound, Mouth,
Skin

Anaerobic Culture Swab (with Gel)



Ordered test(s):

Anaerobic Culture

BD Affirm VPIII



Ordered test(s):

*Vaginitis Panel/DNA Probe
(Trichomonas, Gardnerella,
Candida Species)*

Liquid Based Pap Tests

Ordered Test(s):

*Pap
HPV
Chlam/GC*

*Urine Cytology
Trichomonas*



ThinPrep



Brush



Broom



Paddle

Chlam / GC AND Trichomonas



Unisex Genital Swab

OR



Unisex Urine Collection



ThinPrep

Biopsy Specimens



10% Formalin

*Ordered test(s):
Biopsy*